



HANDLE GUIDE

** Always measure with the handle out of the furniture **



Large oval handle -
Available in Black and Brown plastic, Black and Chrome effect metal - 6 mm Barrel attachment only



Small oval handle -
Available in Black and Brown plastic - 3.5 mm Barrel or 6 mm Barrel attachment

INSTRUCTIONS TO MEASURE & IDENTIFY HANDLE:

1. Completely remove handle via screws located top and bottom
2. Completely remove cable from back of handle
3. Measure A
4. Measure B
5. Measure C
6. Measure cable if required
7. Order from website

** All measurements are a guide only **

Call our Tech Support Line on 0488 603 876

PO Box 1577, North Lakes Qld 4509 | info@sofarehab.com.au | www.sofarehab.com.au | ABN 73 620 447 565