

## HANDLE GUIDE

\*\* Always measure with the handle out of the furniture \*\*





## Large oval handle -Available in Black and Brown plastic, Black and Chrome effect metal - 6 mm Barrel

attachment only



## Small oval handle -

Available in Black and Brown plastic - 3.5 mm Barrel or 6 mm Barrel attachment

## **INSTRUCTIONS TO MEASURE & IDENTIFY HANDLE:**

- 1. Completely remove handle via screws located top and bottom
- 2. Completely remove cable from back of handle
- 3. Measure A
- 4. Measure B
- 5. Measure C
- 6. Measure cable if required
- 7. Order from website

\*\* All measurements are a guide only \*\*

Call our Customer Support Line on 07 2111 4088

130 mm

PO Box 1577, North Lakes Qld 4509 | info@sofarehab.com.au | www.sofarehab.com | ABN 73 620 447 565

